

**MANUFACTURED
IN THE UNITED STATES**

US Patent Pending

**Using only high quality internal
and external materials.**

**Exclusively designed for
Comfort and Purpose.**

**An emergency device that works
outside or inside the home.**

Security practically anywhere.

The back of this device and the charging cable both contain magnets. If you have a health device susceptible to magnetic fields, such as a pacemaker or defibrillator, please check with your physician prior to using this device.

We highly recommend that you do NOT wear this device on a pendent or keep it in a shirt pocket if you have a pacemaker or a defibrillator.



MedEquip Alert™

A Personal Emergency
Response System



www.MedEquipAlert.com

1-844-HELP-808
(4357)



How The MedEquip Alert™ Works:

Easy setup; one recessed emergency SOS call button.

Pushing this button and then releasing it begins a call to a monitoring center. Speak directly to the monitoring center through your **MedEquip Alert™**. The unit will power off when the call is ended.

Audio – Features an audio high quality two-way speaker, you talk directly to the Monitoring center for clear sound to be heard by you and the center’s operator.

The **MedEquip Alert™** will include voice prompts that inform you the device is charging, the battery is low, is calling the emergency center, or a call has been ended.

Charging & Battery Life:

Uses a rechargeable battery and magnetic charging cable. Can be charged in car or plugged into cell phone charger or USB port.

Includes a USB charging cord for easy recharging.

Stand by power will allow the unit to operate on a single charge for 48 hours.

Supports 30 minutes of talk time from a full charge.

Additional Features:

Water proof to the industry standard. It can be worn in the shower, bath, etc.

Roams seamlessly throughout the USA.

Operates from 0 degrees to 125 degrees F.

Optional LOCK BOX Available:

MedEquip Supply will supply a **lock box for a key** to the home entrance doorway. This allows the emergency team to enter without causing extensive damage when arriving at your home.

Product Information:

Does not require a landline. Suitable for all emergencies, accidents, illnesses or dangers.

Wear around the neck or clipped to your belt or pocketbook for 24-hour, 7-day emergency protection.

The MedEquip Alert™ will periodically signal you with information such as battery level and cellular signal strength. The user will also be prompted to recharge the device whenever the battery level falls below 20%.

The MedEquip Alert™ is designed that, if necessary, the firmware in the handset can be upgraded over the wireless cellular network, giving you the most up-to-date means of communication.



Cellular network is AT&T LTE

GPS is included; emergency services will know your location

Included:

Lanyard
Belt Clip
USB Magnetic Charging Cable and Socket
Operating Instructions

Contact Information:

Local Office:
516- 477-6523

National Contact:
1-844-HELP-808
(4357)

Fax:
516-484-1982

www.medequipalert.com