

**Med Equip Alert**  
**Designed for**  
**America's Veteran**



**GPS ACTIVE**  
**SAFETY AND FREEDOM**

**NO MONTHLY FEES**

**NO SUBSCRIPTIONS FEES**

**FEAR OF FALLING?**

**FALL DETECTION**  
**AVAILABLE UPON**  
**REQUEST**

**MANUFACTURED**  
**IN THE UNITED STATES**  
Joint Commission Accredited

**GSA APPROVED**  
**GS-35F-202GA**

**ITEM NUMBER-MA**

**OPERATIONAL COAST**  
**TO COAST**

**WATER RESISTANT**  
**NO LAND LINE NEEDED**  
**CHARGEABLE IN YOUR CAR**  
**SIMPLE TO CHARGE AS YOUR**  
**CELL PHONE, NO CRADLE**  
**REQUIRED**

**Have peace of mind**  
**be safe and secure**  
**7 days a week 24 hours a day**



**Medequipalert.com**

**BRING THIS FLYER**  
**TO YOUR VA**  
**PRIMARY CARE**  
**DOCTOR OR VA**  
**CLINICIAN AT THE**  
**VA HOSPITAL YOU**  
**VISIT AND HAVE**  
**THEM SUBMIT A**  
**CONSULT TO**  
**PROSTHETICS FOR**  
**The Med Equip Alert™**  
**INCLUDING THE**  
**ITEM NUMBER-MA**  
**AND GS-35F-202GA**

**Please include with Fall**  
**Detection, if needed**

**Contact Information**

**Main Office: 844-HELP-808**  
**Marketing: 516-477-6523**  
**Fax: 516-484-1982**

**CUSTOMER SERVICE: To activate**  
**your ALERT: 1-800 327-7010**

## OPTIONAL LOCK BOX AVAILABLE

Med Equip Supply can supply **A LOCK BOX** for a **KEY** to the home entrance doorway. This allows the emergency team to enter without causing possible extensive damage when arriving at your home.

### Charging and Battery Life:

The **Med Equip Alert™** uses a rechargeable battery and includes a **USB Charging cord for easy recharging that can also be used in your car with a USB port.**

Connects to a Monitoring Center, all cellular and GPS antennas are internal.

Easy setup, ONE recessed Emergency

SOS call button.

**Pushing the button for 2 seconds and then releasing it begins a call to a Monitoring center. To cancel the call press and hold for 3 seconds.**

**The MA will include voice prompts that inform you the device is charging, or a call is being placed, or a call has ended.**

**Audio** –features high quality two-way speaker phone, you talk directly to the Monitoring center.

**Water Resistant** to the industry standard.


## FREEDOM TO GO OUTSIDE WHEN YOU WANT

The **Med Equip Alert™** allows the user freedom to leave the house and enjoy independent living 24 hours, 7 days a week security and protection, putting you in touch with an emergency operator that has your necessary information.

## **FEATURE PLATFORM**

**STAT-VA-DATA** connects the veteran's medical records to the Med Equip Alert by scanning the QR code on the back of the alert. All of the information is shared with first responders, doctors and any other medical clinicians. This will allow them to start emergency treatment without delay based on this information.

**TRY IT OUT  
SCAN THE QR CODE!**



**Select yes, then select any of the topics-this could be your information for them to see.**

**Roams seamlessly** throughout the USA Suitable for all emergencies, accidents, illness or danger.

## REMOTE MANAGEMENT SIGNALING

The **Med Equip Alert™** is designed, that the Firmware in the Handset can be upgraded over the wireless cellular network giving you longevity of the alert and the most up to date means of communication.

## INCLUDED WITH THE Med Equip Alert™

Lanyard  
Belt Clip

USB Magnetic Charging Cable and Socket

The Med Equip Alert uses cellular and or digital network to send and receive information if this network is not receiving the information for any reason the data will not be received by the monitoring service center, reasons are numerous such as the equipment is out of range because of location, building blockage or weather conditions. The device may not always work because of the above reasons. Med Equip is not responsible if the device cannot connect to the proper network.

The back of this device and the charging cable both contain magnets. Magnets may cause permanent damage to health devices susceptible to magnetic fields, credit cards, computer hard drives, watches, TVs, data storage media and other electronic devices. If you have a health device susceptible to magnetic fields, such as a pacemaker or defibrillator, please check with your physician prior to using this device.

We highly recommend that you do NOT wear this device on a pendant or keep it in a shirt pocket if you have a pacemaker or a defibrillator. The warranty is void if the product has been damaged or tampered with or any such parts have been opened. Limited warranty covers defects in materials and workmanship in your Med Equip Alert for one year from date of purchase. If user does not show intended use for a period of six months service will be suspended. User must call customer service to reactivate the device.

Patent US D 871,359S